

AESTHETICS



Ammad Surgical

Vision to precision



MERCEDES

Three openings in a circumferential pattern



ACCELERATOR

Three openings in triangle pattern. One center hole and two lateral holes



CANE

Two elongated holes



CANDY CANE

Three elongated holes



TWO PORT RADIAL

Two lateral holes



Harvest four holes

Two lateral holes and two holes on opposite side



STANDARD ONE HOLE

One hole near the tip



STANDARD TWO HOLES

Two holes near the tip in linear pattern



STANDARD THREE HOLES

Three holes in linear Pattern



SATTLER

Twenty four holes. Eight on each side



SPATULA

One hole near the tip



TOLEDO V

One hole

INFILTRATION CANNULA



Eight holes



Ten holes



Micro facial




Luer lock II



Luer lock III



Luer lock I



Harvest one hole



Half monty twenty seven holes



Full monty fifty four holes

Fat Transfer



Aluminum



Stainless Steel

INFILTRATION

- Four Hols
- Eight Holes
- Ten Holes
- Micro Facial
- Micro Cannulas

ASPIRATION

- Mercedes
- Accelerator
- Cane
- Candy Cane
- Two Port Radial
- Harvest Four Holes
- Standrad One Hole
- Standrad Two Holes
- Standrad Three Holes
- Sattler
- Spatula
- Toledo V



Disposable Cannulas are packed in double peel away pouches Srerlised (EO) as shown


DISPOSABLE LUER LOCK

Luerlock hub




DISPOSABLE

Single piece cannula



Note: For custom made cannulas you can send us drawings / specifications at:
info@ammadsurgical.com

CANNULA SIZES




CANNULA TIP STYLES



CURVE PREFERENCE

- a. Convex Curve
- b. Straight
- c. Concave Curve
- d. "S" TIP





Ammad Surgical
Vision to precision

CANNULA HANDLES AND PACKING

COLOURS AVAILABLE




LUER/LUER

Female Luer on top for
Cannulas connections &
male luer connection for
syringe

ALUMINUM HANDLE

ANODISED & Knurled



TUBE PACKING

Head Office

41-B Commercial Area,
Cavalry Ground, Cantt.,
Lahore, Pakistan
Tele: +92-(0)42-66-100-70,66-807-67
FAX: +92-(0)-66-679-05
info@ammadsurgical.com

Factory

Defence Road,
Noorabad,
Sialkot, Pakistan
Tele: +92-(0)52-32-566-59

