

DIGNOSTIC SET

◀ ENT DIGNOSTIC SET ▶

Details:

- Boyonet locking mounts for each instruments head to allow secure connection to battery handle
- Rheostatic ON / OFF Switch
- Otoscope with 3 Ear specula & Nasal speculam with mannual screw opening
- Ophthalmoscope head with lens wheel featuring 25 corective lenses
- Chrome plated bent arm larynx with light
- 2 Mirrors 4mm & 5mm
- Disposable tongue depressor with holder.



PLASTIC PACKING



LEATHER POUCH

Fiber optic
McCoy



Fiber optic laryngoscope

Stainless steel



Disposable blade with fiber optic light

DIGNOSTIC SET



◀ ENT DIAGNOSTICS SET (SCREW)

Details:

Screw locking mounts for each instrument head to allow secure connection to the battery handle "C" Cell powered. Stain Finish chrome plated batter handle Rheostatic on / off switch Total attachments 11

1. Otoscope:

i. Chrome plated otoscope head with removable 4x magnifying lens and inufflator fitting ii. Autoclavble polypropylene ear specula in 2.5mm, 3.5 mm and 4.5 mm iii.

Nasal speculum adaptor with manual screw opening for instrument connected to the otoscope head

2. Ophthalmoscope:

i. Ophthalmoscope head with lens wheel featuring 25 corrective lens from -25 to +40 diopters

3. Bent Arm Illuminator i. Chrome plated bent arm larynx with a bulb ii. One 4.0 mm and one 5.0 mm mirror iii.

One tongue depressor adaptor

Stainless Steel iv. One tongue depressor



◀ ENT DIAGNOSTICS SET (LOCK)

Details:

Boyonet locking mounts for each instrument head to allow secure connection to the battery handle "C" Cell powered. Stain Finish chrome plated batter handle Rheostatic on / off switch Total attachments 11.

Otoscope: Chrome plated otoscope head with removable 4x magnifying lens and inufflator fitting Autoclavble polypropylene ear specula in 2.5mm, 3.5 mm and 4.5 mm Nasal speculum adaptor with manual screw opening for instrument connected to the otoscope head

2. Ophthalmoscope: Ophthalmoscope head with lens wheel featuring 25 corrective lens from -25 to +40 diopters

3. Bent Arm Illuminator Chrome plated bent arm larynx with a bulb

One 4.0 mm and one 5.0 mm mirror

One tongue depressor adaptor

Stainless steel One tongue depressor

www.ammadsurgical.com

DIAGNOSTIC SET



OTOSCOPE FIBER OPTIC (Pocket) ▶

Details:

- 2.5 V Halogen Light
 - Fiber Optic for Obstruction Free Light Sealed System with Insufflator Fitting for Pneumatic Otoscope
 - Otoscope Head Made from Polymid
 - 7 Disposable Polypropylene Ear Specula
- Colours:** Black, Blue, Pink, Green.



◀ OTOSCOPE

Details:

- Boyonet locking mounts for each instruments head to allow secure connection to battery handle
- Rheostatic ON / OFF Switch
- Otoscope with 3 Ear specula



◀ OPHTHALMOSCOPE

- Boyonet locking mounts for each instruments head to allow secure connection to battery handle
- Rheostatic ON / OFF Switch
- Ophthalmoscope head with lens wheel featuring 25 corrective lenses