

# DIGNOSTIC SET

## ◀ ENT DIGNOSTIC SET ▶

Details:

- Boyonet locking mounts for each instruments head to allow secure connection to battery handle
- Rheostatic ON / OFF Switch
- Otoscope with 3 Ear specula & Nasal speculam with manual screw opening
- Ophthalmoscope head with lens wheel featuring 25 corective lenses
- Chrome plated bent arm larynx with light
- 2 Mirrors 4mm & 5mm
- Disposable tongue depressor with holder.



PLASTIC PACKING



LEATHER POUCH

Fiber optic  
McCoy



Fiber optic laryngoscope



Stainless steel



Disposable blade with fiber optic light

# DIGNOSTIC SET



## ◀ ENT DIAGNOSTICS SET (SCREW)

Details:

Screw locking mounts for each instrument head to allow secure connection to the battery handle "C" Cell powered. Stain Finish chrome plated batter handle Rheostatic on / off switch Total attachments 11

1. Otoscope:
  - i. Chrome plated otoscope head with removable 4x magnifying lens and inufflator fitting
  - ii. Autoclavble polypropylene ear specula in 2.5mm, 3.5 mm and 4.5 mm
  - iii. Nasal speculum adaptor with manual screw opening for instrument connected to the otoscope head
2. Ophthalmoscope:
  - i. Ophthalmoscope head with lens wheel featuring 25 corrective lens from -25 to +40 diopters
3. Bent Arm Illuminator
  - i. Chrome plated bent arm larynx with a bulb
  - ii. One 4.0 mm and one 5.0 mm mirror
  - iii. One tongue depressor adaptor
- iv. One tongue depressor



## ◀ ENT DIAGNOSTICS SET (LOCK)

Details:

Boyonet locking mounts for each instrument head to allow secure connection to the battery handle "C" Cell powered. Stain Finish chrome plated batter handle Rheostatic on / off switch Total attachments 11.

1. Otoscope: Chrome plated otoscope head with removable 4x magnifying lens and inufflator fitting Autoclavble polypropylene ear specula in 2.5mm, 3.5 mm and 4.5 mm Nasal speculum adaptor with manual screw opening for instrument connected to the otoscope head
2. Ophthalmoscope: Ophthalmoscope head with lens wheel featuring 25 corrective lens from -25 to +40 diopters
3. Bent Arm Illuminator Chrome plated bent arm larynx with a bulb
- One 4.0 mm and one 5.0 mm mirror
- One tongue depressor adaptor
- Stainless steel One tongue depressor

# DIGNOSTIC SET



## OTOSCOPE FIBER OPTIC (Pocket) ▶

Details:

- 2.5 V Halogen Light
  - Fiber Optic for Obstruction Free Light Sealed System with Insufflator Fitting for Pneumatic Otoscope
  - Otoscope Head Made from Polymid
  - 7 Disposable Polypropyleni Ear Spcula
- Colours:** Black. Blue. Pink. Green.



## ◀ OTOSCOPE

Details:

- Boyonet locking mounts for each instruments head to allow secure connection to battery handle
- Rheostatic ON / OFF Switch
- Otoscope with 3 Ear specula



## ◀ OPHTHALMOSCOPE

- Boyonet locking mounts for each instruments head to allow secure connection to battery handle
- Rheostatic ON / OFF Switch
- Ophthalmoscope head with lens wheel featuring 25 corective lenses