

HOLLOW WARE

To increase the working life of our products we use Rotary Vibrators for the De Burring and polishing of the instruments. Using this technology we are able to polish our instruments in areas where hand polishing is not possible giving us an edge over our competitors.



FORCEP JAR



STERILIZING DRUM

Size: 120 x 100mm
150 x 150mm
190 x 145mm
190 x 240mm
240 x 145mm
240 x 240mm
390 x 390mm



KIDNEY TRAY

Size: 6"
8"
10"



INSTRUMENT BOX

Size: 180 x 80 x 40mm
250 x 120 x 60mm
400 x 200 x 75mm
500 x 200 x 120mm